CERTIFICATE OF ANALYSIS

Complement Technology, Inc. 4801 Troup Hwy, Suite 701 Tyler, Texas 75703, USA	Product: Hemolysin Cat#: Hemolysin Lot#: Expiration Date:
Description: Hemolysin	
Specifications	
FILL VOLUME	2.1 mL
COMPOSITION	Rabbit anti-sheep erythrocyte antiserum in 50 % Glycerol.
PRESERVATIVE	None
STORAGE CONDITIONS:	Hemolysin should be stored at 4 °C during usage. Hemolysin may be stored for longer term at ≤ -20°C
Avoid Repo	eated Freeze/Thaw.
	RESEARCH USE ONLY JMAN OR DRUG USE
Signature of Analyst	 Date

Phone: 1-903-581-8284 FAX: 1-903-581-0491 Email: contactCTI@ComplementTech.com

Web Site: www.ComplementTech.com